



# ***Digital Preoperative Planning for Long-Bone Fractures***

Esmitt Ramírez

Ernesto Coto

[esmitt.ramirez@ciens.ucv.ve](mailto:esmitt.ramirez@ciens.ucv.ve)

[ernesto.coto@ciens.ucv.ve](mailto:ernesto.coto@ciens.ucv.ve)

Universidad Central de Venezuela




# Introduction

- Nowadays, digital images are a fundamental base to medical diagnosis
- CAD (Computer Aided Diagnosis) systems provide support in preoperative surgery planning
- CAOS – Computer Aided Orthopedic Surgery




# Preoperative Planning

- It helps estimating the final result of the surgery
- It consists of a trace of the fracture, the steps to follow during surgery and the implant to be used
- Often the tracing is done manually in paper



# Preoperative Planning



## Pre-Op Plan

1. Supine/Arm table / Tourniquet No Exsanguination
2. "Henry's" (Anterior) Approach to Radius
3. Expose fx Ø stripping "Specials": Dental pick
4. Reduce → Interfragmentary screw
5. Neutralization plate - 7 hole LCD plate (leave central hole open)
6. Ulnar - approach to ulna: ECU/FCU interval
7. Reduce Make sure no comminution on interosseous membrane!
8. Compression Plate
  - Screw #1 - Hole 6 (proximal) in neutral
  - Screw #2 - Hole 2 (distal) in compression
  - Fill holes 1/3/5/7 neutral
9. ✓ Fluoro + full ROM
10. Close subcut d skin (Drain if needed)
11. No splint



# CAOS system

- Manual procedure (pen-paper) consumes long time
- Manual procedure is error-prone
- Instead, CAOS systems are used in many healthcare centers
- CAOS systems make **digital planning** possible

**Then, our proposal is**


**Eliminate the pen-paper method**

**Build a low-cost CAOS solution**

**Special hardware must not be required**



# Proposed Scheme






# Image Acquisition



## Prerequisites

- Controlled environment (lighting, distance to shoot)
- Resolution of at least 800x600 pixels (0.5 megapixels)
- Focus in the area of interest (i.e. the fracture)

# Image Calibration




Distance between holes is configurable

Standard ISO 838 establishes this distance is  $12 \pm 1$  mm



# Image Enhancement

- Necessary to obtain a better image contrast
- Manual window/level adjustment



# Segmentation and Puzzling

- **Segmentation:** Extract from the image all bone segments separated by the fracture
- **Puzzling:** Place these bone segments in their anatomically correct positions using translations or rotations

## Segmentation



Manual

Semi-automatic

In manual mode, the surgeon has full control over the segmentation

In semi-automatic mode, a Canny border detection algorithm is used to find the border of each bone segment

# Segmentation and Puzzling



A fracture would produce many bone fragments


The measurement tool is very useful to calculate the length of the bone, the length of the implant, etc.

For a later reviewing, different annotations should be written

Not all segments fit exactly in their correct anatomical position, a very close approximation is enough for planning

# Implant Placement

- When the puzzling process is completed, the surgeon decides whether to put an implant on the patient or not
- A library of implants is provided, including: plates, screws and pins
- MySql® is used to store information about implants
- Each implant is represented using a STL file which stores the 3D triangulated surface of the implant
- A complete library of AO<sup>2</sup> fractures is also available





# Implant Deformation

- In some cases, surgeons might need to bend an implant
- The implant deformation process is useful to plan the bending operation over the implant during the surgery
- At the surgery room, the surgeon applies a deformation over the implant using a specialized groove plier
- Our system allows the digital deformation of the implant prior to the surgery



# Tracing Paper and Report Generation

- Finally, the system can print out the tracing paper and a report that includes: annotations, measurements, fragments of bones, etc.

**At this point, the  
preoperative planning  
is complete!**


HOSPITAL CLINICO DE CARACAS

DEPARTAMENTO DE RADIOLOGIA

DR. CARLOS SANCHEZ

PACIENTE: LUIS PERRAZO EDAD: 45 PESO: 78 FECHA: 01/11/2009

ENTRADA: PACIENTE RECIBIDO POR FRACTURA TIPO 32-A3



TACTICA A EMPLEAR:

QUITAR EL TORNIQUETE

LIMPIAR LA ZONA

APLICAR CORTE TIPO 4  
BAJO ANESTESIA LOCAL

SELECCIONAR EL  
IMPLANTE ESCOGIDO

NO APLICAR  
DOBLAMIENTO

RECONOCER LOS  
LIGAMENTOS COMO AREA  
DE SEGURIDAD

COLOCAR EL IMPLANTE SIN  
ARRAIGAMIENTO PARCIAL  
SINO TOTAL

REQUIERE DE IMPLANTE EXTERNO. UNA PLACA DEL TIPO DCP ANATOMICA  
MARCA XXX COLOCADA A 5MM DEL SEGMENTO DIAFISIAL



# Results and Conclusions

- We installed our system in a standard PC without specialized hardware at the Radiology Department of the HUC (University Hospital of Caracas) in Venezuela
- Two experts were instructed on how to use the system in a short time
- Without the system, the surgeons create a preoperative planning manually in 8 to 25 minutes approximately
- With the system the planning time varied from 2 to 5 minutes
- That represents a considerable reduction of the planning time and positive feedback about results of the system



# Future works

- Implementing a better border detection algorithm, e.g. algorithms based on active contours or deformable models
- Extends our system to support others object of reference
- Take advantages of current technology to speed up the annotation process, e.g. using speech recognition tools, medical tablet PC's, etc.

# PLANIFICACIÓN PREOPERATORIA

TRAUMATOLOGÍA Y ORTOPEDIA  
VERSIÓN 1.0 - 2009

ESMITT RAMIREZ, CARLOS SANCHEZ



**Realizar  
Planificación**



**Clasificación  
AO**



**Biblioteca  
de Implantes**



**Generar  
Reportes**



**Base de  
Datos**



**Configuración**



**Modificar**



**Salir**



# Questions



<http://ccg.ciens.ucv.ve>

{esmitt.ramirez | ernesto.coto}@ciens.ucv.ve